

Targeting Systems Introduces : New Products

Innovative tools for studying gene expression

New Products : psiScreen system- Gaussia luciferase-based expression vectors and assay kits for rapid and quantitative evaluation of gene silencing.

The psiScreen Vectors are designed to provide a rapid, quantitative approach for evaluation and optimization RNA interference (RNAi). These vectors enable monitoring of changes in expression of a target gene fused to a novel luciferase reporter gene. In both vectors, Gaussia luciferase is used as a primary reporter gene, and the target gene of interest can be cloned into multiple cloning sites located downstream of the translational stop codon of the luciferase gene. Transfection of cells with the psiScreen vectors results in the production of an mRNA in which the mRNA encoding Gaussia luciferase is expressed as a fusion with the mRNA encoding the target gene. Initiation of the gene silencing occurs when co-transfection of siRNA towards the target gene of interest results in cleavage and subsequent degradation of fusion mRNA. Measurement of decreased Gaussia luciferase activity serves as an indicator of RNA interference. Since the Target gene is sub cloned after the stop codon of Gaussia luciferase, the activity of Gaussia luciferase is unaffected by the fusion partner.

Components of the psiScreen system

Gaussia Luciferase-based expression vector and reagents for siRNA screening

Catalog #	Capacity	Price
#S001	25 µg psiScreen, plus GAR-1 reagent (100 assays)	\$420
#C001	pCMV-Fluc (25 µg)	\$200

Gaussia Luciferase RNAi Assay Kit

Product cord/code	Capacity	Price
GAR-1	1000 assays	\$350
GAR-2	1000 assays	\$350

Advantages

An easy method to screen a wide variety of siRNAs for gene silencing :

The effectiveness of different siRNAs to silence the target gene of interest is evaluated quantitatively simply by measuring the luciferase activity.

Simplicity :

Gaussia luciferase is secreted into the media. It is therefore necessary to only assay cell supernatants for luciferase activity without the need for lysing the cells. Considerable time is saved since time course experiments can be performed using the same group of transfected cells without lysing at each time point.


Sensitivity :

Gaussia luciferase, a thermostable enzyme, is 1000 times brighter than Renilla and firefly Luciferase thus increasing sensitivity of the assay.

Convenience of an easy-to use assay kit :


Using the GAR reagents, a simple one-step reporter assay is possible.

The GAR-2 reagent has stabilizing components which stabilize the bioluminescent signal enabling the system to be useful for high throughput screening applications



The psiScreen vector

The psiScreen vector has a multiple cloning site after the stop codon of the Gaussia luciferase gene



Features of the psiScreen-Basic vector :

- pCMV-GLUC-1 (5764 bp)
- CMV promoter bases: 209-863
- Gaussia luciferase gene: 907-1497
- T7 promoter bases: 864-882
- Polylinker bases: 889-907
- SP6 promoter: 1513-1530
- Synthetic polyadenylation site: 1497-1541
- SV40 promoter bases: 2082-2417
- SV40 origin of replication: bases 2196-2281
- Neomycin ORF : bases 2453-3247
- SV40 PolyA: bases 3302-3674
- ColE1 origin: bases 3934-4607

The sequence of the psiScreen vector can be found at our website www.targetingsystems.com

And also on the last page of this brochure




Figure 2 : This figure shows dose response of siRNA for gene silencing of the target gene (human tumor suppressor p53 gene). The data has been normalized for transfection efficiency using a firefly luciferase expression vector as a control plasmid. Unrelated siRNA was shown to be ineffective in gene silencing (data not shown).

Evaluation of the siScreen system for screening siRNAs against the human tumor suppressor p53 gene :

Different siRNAs against human p53 gene were synthesized and evaluated for their ability to silence the p53 gene by co-transfection with a siScreen vector in which the p53 coding sequence was subcloned after the stop codon. The sequences and scoring of the siRNA duplexes used for evaluation are shown below. As expected siRNA sequences with a high score showed very effective gene silencing of the target genes. The data p-lotted shows results after normalization for transfection efficiency using a firefly luciferase expression vector pCMV-Fluc as a denominator plasmid.




Figure 3 : This figure shows the effect effectiveness of different siRNAs in silencing the p53 gene. HEK'293 cells were cotransfected with the siScreenp53 expression vector and siRNA against either the p53 gene oligo 2 , and siRNA mix form NEB New England Biolabs, siRNA against gaussia luciferase oligo 1 or unrelated siRNA green bar . The cell supernatants were assayed for gaussia luciferase activity 24 hrs after transfection. The results have been normalized for transfection efficiency using a firefly luciferase expression vector HIV'Luc as a denominator plasmid. Transfection

complexes were prepared by mixing 6 µg of pCMV⁺Gluc⁺p53 vector with 3 µg pHIVFluc and 100 pmols siRNA with 12 µl of the Targefect F2 reagent in 0.5 ml of high glucose DMEM. The contents of the tube were mixed well after each addition and incubated at 37° C for 30 min to form the transfection complexes. 65 µl of the transfection complexes were added to 0.5 ml of the supernatant cell culture media containing 5% serum to cells in 24'well dishes. The dishes were gently swirled and then incubated at 37° C overnight and assayed for luciferase activity at 24 hrs post-transfection.

Sequence of siScreen vector with MCS site

```

1 gacggatcgg gagatctccc gatcccctat ggtcgactct cagtacaatc tgctctgatg
61 ccgcatagtt aagccagtat ctgctccctg cttgtgtgtt ggaggtcgct gagtagtgcg
121 cgagcaaaat taaagctaca acaaggcaag gcttgaccga caattgcatg aagaatctgc
181 ttagggtag gcgtttgcg ctgcttcgag atgtacgggc cagatatacg cgttgacatt
241 gattattgac tagtattaa tagtaataca ttacggggtc attagtcat agcccatata
301 tggagttccg cgttacataa cttacggtaa atggcccgc tggctgaccg cccaacgacc
361 cccgccatt gacgcaata atgacgtatg ttcccatagt aacgccaata gggactttcc
421 attgacgta atgggtggac tatttacggt aaactgcca cttggcagta catcaagtgt
481 atcatatgcc aagtacgcc cctattgacg tcaatgacgg taaatggccc gcctggcatt
541 atgccagta catgacctta tgggacttcc ctacttgga gtacatctac gtattagtca
601 tcgctattac catggtgatg cggtttggc agtacatcaa tgggctgga tagcggtttg
661 actcacgggg atttccaagt ctccaccca ttgacgtaa tgggagttg tttggcacc
721 aaaatcaacg ggactttcca aaatgtcgt acaactccgc ccattgacg caaatgggag
781 gtaggcgtgt acgggtggag gtctatataa gcagagctct ctggctaact agagaacca
841 ctgcttactg gcttatcga ataatatcga ctactatag ggagaccaa gcttggtacc
901 gagctcggat ccagccacca tgggagtaa agttctgtt gccctgatct gcatcgctgt
961 ggccgaggcc aagcccaccg agaacaacga agactcaac atcgtggccg tggccagcaa
1021 cttcgcgacc acggatctcg atgctgaccg cgggaagttg cccggcaaga agctgccgct
1081 ggaggtgctc aaagagatgg aagccaatgc ccggaaagct ggctgcacca ggggctgtct
1141 gatctgctg tcccatca agtgacgcc caagatgaag aagttcatcc caggacgctg
1201 ccacacctac gaagcgaca aagagtcgc acagggcggc ataggcgagg cgatcgctga
1261 cattcctgag attcctgggt tcaaggactt ggagccatg gagcagtca tcgcacaggt
1321 cgatctgtgt gtggactgca caactggctg cctcaaaggg cttgccaacg tgcagtgttc
1381 tgacctgctc aagaagtggc tgccgcaacg ctgtgcgacc ttgcccagca agatccaggg
1441 ccaggtggac aagatcaagg gggccggtgg tgactaagcg gccgctcgag agatcttggg
1501 attctgcaga tatctctaga aataattctt actgtcatgc caagtaagat gcttttctgt
1561 gctgcaatag caggcatgct ggggatgctg tgggctctat ggcttctgag gcggaaagaa
1621 ccagctgggg ctctaggggg tatcccacg cgccctgtag cggcgcatta agcgcggcgg
1681 gtgtggtggt tacgcgacg gtgaccgcta cacttgccag cgccctagcg cccgctcctt
1741 tcgcttctt ccttcttct ctgccacgt tcgcccgtt tcccgtcaa gctctaaatc
1801 ggggcatccc tttaggggtc cgatttagtg cttacggca cctcgacccc aaaaaactg
1861 attaggtgta tggttcacgt agtgggcat cgccctgata gacggttttt cgcccttga
1921 cggtggagtc cacgttctt aatagtggac tctgttcca aactggaaca aactcaacc
1981 ctatctcggg ctattctttt gattataag ggatttggg gatttcggcc tattggttaa
  
```

2041 aaaatgagct gattaacaa aaattaacg cgaattaatt ctgtggaatg tgtgtcagtt
2101 aggggtgtgga aagtccccag gctccccagg caggcagaag tatgcaaagc atgcatctca
2161 attagttagc aaccaggtgt ggaaagtccc caggctcccc agcaggcaga agtatgcaaa
2221 gcatgcatct caattagta gcaacatag tcccgccct aactccgcc atcccccccc
2281 taactccgcc cagttccgcc cattctccgc cccatggctg actaatttt ttatttatg
2341 cagaggccga ggccgctct gcctctgagc tattccagaa gtagtgagga ggctttttg
2401 gaggcctagg cttttgcaa aagctcccg gagctgtat atccatttc ggatctgatc
2461 aagagacagg atgaggatcg ttccgatga tgaacaaga tggattgcac gcaggttctc
2521 cggccgctg ggtggagagg ctattcggt atgactggg acaacagaca atcggctgct
2581 ctgatccgc cgtgtcccg ctgtagcgc aggggcccc ggttctttt gtcaagaccg
2641 acctgtccg tgccctgaat gaactgcagg acgaggcagc gcggctatc tggctggcca
2701 cgacggcgct tcctgcgca gctgtgctg acgtgtcac tgaagcggga agggactggc
2761 tgctattgg cgaagtccg gggcaggatc tcctgtatc tcacctgtc cctgccgaga
2821 aagtatccat catggctgat gcaatggcg ggctgcatac gcttgatccg gctacctgcc
2881 cattcgacca ccaagcгаа catgcatcg agcagcagc tactcggatg gaagccggtc
2941 ttgtgatca gtagatctg gacgaagagc atcaggggct cgcgccagcc gaactgttcg
3001 ccaggctcaa ggcgcgatg cccgacggcg aggatctct cgtgacctat ggcgatgcct
3061 gcttgccgaa tatcatggtg gaaaaaggcc gcttttctgg atcatcgac tgggcccggc
3121 tgggtgtggc ggaccgctat caggacatag cgttggctac ccgtgatatt gctgaagagc
3181 ttggcggcga atgggctgac cgcttctcg tgctttacgg taccgccgct cccgattcgc
3241 agcgcacgc ctctatcgc ctcttgacg agttctctg agcgggactc tggggttcga
3301 aatgaccgac caagcgacgc ccaacctgcc atcacgagat ttcgattcca ccgccgcctt
3361 ctatgaaagg ttgggctcg gaatgcttt ccgggacgcc ggctggatga tcctccagcg
3421 cggggatctc atgctggagt tcttcgcca cccaacttg ttattgtag ctataatgg
3481 ttacaaataa agcaatagca tcacaaatt cacaaataa gcatttttt cactgcattc
3541 tagttgtgt ttgtcctaac tcatcaatg atcttatcat gtctgtatc cgctgcacctc
3601 tagctagagc ttggcgaat catggtcata gctgttctc gtgtgaaatt gttatccgct
3661 cacattcca cacacatac gagccggaag cataaagtgt aaagcctggg tgcctaatg
3721 agtgagctaa ctacattaa ttgcgtgag ctactgccc gctttccagt cgggaaacct
3781 gtcgtgccag ctgcattaat gaatcggcca acgcgcgggg agaggcggtt tgcgtattgg
3841 gcgctctcc gcttctcgc tcaactgact gctgcgctc gtcgttcggc tgcggcgagc
3901 ggtatcagct cactcaaagg cggtaatagc gttatccaca gaatcagggg ataacgcagg
3961 aaagaacatg tgagcaaaag gccagcaaaa ggccaggaac cgtaaaaagg ccgctgtgct
4021 ggcgttttc cataggctcc gccccctga cgagcatcac aaaaatcgac gctcaagtca
4081 gaggtggcga aaccgacag gactataaag ataccaggcg ttccccctg gaagctccct
4141 cgtgctctc cctgtccga cctgcccgt taccggatac ctgtccgct ttctccctc
4201 ggggaagcgtg gcgctttctc aatgctcacg ctgtaggtat ctgattcgg ttaggtcgt
4261 tcgctccaag ctgggctgtg tgcacgaacc cccgctcag cccgaccgct gcgccttctc
4321 cggtactat cgtctgagt ccaacctggg aagacacgac ttatgccac tggcagcagc
4381 cactggtaac aggattagca gagcgaggta ttagggcgt gctacagagt tctgaagtg
4441 gtggcctaac tacggctaca ctagaaggac agtatttgg atctgcgctc tgcgtaagcc
4501 agttacctc gaaaaagag ttggtagctc ttgatccggc aaacaaacca ccgctggtg

4561 cgggtggttt tttgttgca agcagcagat tacgcgcaga aaaaaaggat ctcaagaaga
4621 tcctttgatc tttctacgg ggtctgacgc tcagtgaac gaaaactcac gttaagggat
4681 ttggtcatg agattatcaa aaaggatctt cacctagatc cttttaaatt aaaaatgaag
4741 ttttaaatac atctaaagta tatatgagta aactgggtct gacagttacc aatgcttaat
4801 cagtgaggca cctatctcag cgatctgtct atttcttca tccatagttg cctgactccc
4861 cgctgtgtag ataactacga tacgggaggg cttaccatct ggccccagtg ctgcaatgat
4921 accgcgagac ccacgctcac cggctccaga tttatcagca ataaaccagc cagccggaag
4981 ggccgagcgc agaagtggtc ctgcaacttt atccgcctcc atccagtcta ttaattgtg
5041 ccgggaagct agagtaagta gttcggcagt taatagttg cgcaacgttg ttgccattgc
5101 tacaggcatc gtggtgtcac gctcgtcgtt tggatggct tcattcagct ccggttccca
5161 acgatcaagg cgagttacat gatccccat gttgtgcaaa aaagcggta gtccttcgg
5221 tcctccgacg gttgtcagaa gtaagtggc cgcagtgtta tcaactatgg ttatggcagc
5281 actgcataat tctctactg tcatgccatc cgtaagatgc tttctgtga ctggtgagta
5341 ctcaaccaag tcattctgag aatagtgat gcggcgaccg agttgctctt gcccggcgtc
5401 aatacgggat aataccgagc cacatagcag aacttataaa gtgctcatca ttgaaaacg
5461 ttctcgggg cgaaaactct caaggatctt accgctgtg agatccagt cgatgaacc
5521 cactcgtgca cccaactgat cttcagcatc tttacttc accagcgtt ctgggtgagc
5581 aaaaacagga aggcataatg ccgcaaaaaa gggaataagg gcgacacgga aatgttgaat
5641 actcactc ttcttttc aatattattg aagcattat cagggttatt gtctcatgag
5701 cggatacata ttgaaatgta tttagaaaaa taaacaaata ggggttccgc gcacattcc
5761 ccgaaaagtg ccacctgacg tc