

**NEW! UltraBrite See and Assay Kit** (catalog # ULTRA-SA-RedFluc-1)

Visualize with *mVermilion*, quantify with *Luciferase*.

**Kit Components:** LentiGlo plasmids bundled with UltraBrite™ Firefly luciferase assay reagent (FLAR-1)


**How It Works:** A single expression cassette drives co-expression of luciferase and a bright red monomeric RFP (25 Kd, 2X brighter than mCherry) via a self-cleaving 2A peptide. This ensures both signals are expressed efficiently from one promoter.

- **Visual QC:** Red FP fluorescence under standard scopes
- **Quantitative Power:** luciferase readout on any luminometer

**Advantages at a Glance**

- **Visual + Quantitative:** Red FP QC plus luciferase sensitivity
- **Confirms success of transfection** and allows for the non-destructive.
- **Allows real-time monitoring of gene expression**
- **Detect potential issues:** Researchers can visually inspect the cells for overall health and evenness of expression
- **Helpful in high-throughput screens** where many conditions are being tested, the fluorescent reporter can be used as a rapid, early readout to rule out low transfection efficiency or other issues before performing the more labor-intensive luciferase assay.
- **Simple & Fast:** 2-minute lysis; 30–60 min stable signal window
- **Flexible Expression:** Subclonable cassette; alternate FP options, EGFP, tdTomato available
- **Budget-friendly:** Priced competitively with leading suppliers

Catalog#	Product Name	Contents	Price
ULTRA-SA-RedFLuc-1 <i>Special Offer till October 15</i>	UltraBrite Firefly luciferase See and Assay Reagent Kit	LentiGlo plasmid: pLenti-UBC-RedFluc-T2A-mVermilion Firefly luciferase assay reagent (FLAR-1) , 1000 assays	<b>\$800</b> <i>Special Offer till October 15, Regular \$1100</i>
FLAR-1	Firefly luciferase assay reagent	Firefly luciferase assay reagent (1000 assays)	\$400



**FLAR™ Luciferase Assay Reagent — Key Features**

- **Format:** Homogeneous (H) or Non-homogeneous (NH)
- **Sensitivity:** Maximum light output, comparable to Promega Bright-Glo®
- **Signal Stability:** ~30 min (Photinus pyralis) to ~60 min (Luciola italica)
- **Cell Lysis Time:** ~2 minutes
- **Workflow:** Add equal volume of FLAR™ to culture medium (100 µL reagent to 100 µL medium in 96-well) → wait 2 minutes → read

## Protocol

1. Remove 96- or 384-well plates containing mammalian cells from the incubator.
2. Equilibrate cultured cells to room temperature.
3. Add FLAR reagent equal to medium volume (100  $\mu$ L reagent to 100  $\mu$ L medium in 96-well; 30  $\mu$ L in 384-well).
4. Wait 2 minutes for lysis.
5. Measure luminescence in a luminometer.


Signal Half-Life: 30 min (Photinus pyralis) to 1 hr (Luciola italica)

NH format uses 5 $\times$  CLR buffer.

For pre-lysed samples: Mix 5–20  $\mu$ L lysate with 100  $\mu$ L FLAR reagent and read immediately.

**See it before you measure it.** “Experience the UltraBrite™ See+Assay Kit — order now and reduce wasted assays.”

**Other Options: See and Assay options available with our entire panel of single luciferase reporters**-Gaussia luciferase, Renilla luciferase and Cypridina luciferase . Email [info@targetingsystems.net](mailto:info@targetingsystems.net) or call (619) 562 1518 for available options! **Choice of RFP or EGFP co-expressed with a luciferase reporter of your choice from the panel below!**



## Custom Reagents:

We can provide custom formulations to fit your HTS application. Call our tech support team at 1-866-620-4018 or email us [info@targetingsystems.net](mailto:info@targetingsystems.net)